



ELK Biotechnology

NK-3R Rabbit pAb
Catalog Number: ES8524
For research use only.

Overview

Product Name	NK-3R polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthesized peptide derived from NK-3R at AA range; 291-340
Specificity	NK-3R Polyclonal Antibody detects endogenous levels of NK-3R
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	NK-3R
Gene Name	TACR3
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	52kDa
Human Gene ID	6870
Human Swiss-Prot Number	P29371
Alternative Names	Neuromedin-K receptor (NKR) (NK-3 receptor) (NK-3R) (Neurokinin B receptor) (Tachykinin receptor 3)