



ELK Biotechnology

MEK-1 Rabbit pAb
Catalog Number: ES8526
For research use only.

Overview

Product Name	MEK-1 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthetic peptide from human protein at AA range; 330-410
Specificity	The antibody detects endogenous MEK-1 protein
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	MEK-1
Gene Name	MAP2K1 MEK1 PRKMK1
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	43kDa
Human Gene ID	5604
Human Swiss-Prot Number	Q02750
Alternative Names	MAP2K1 MEK1 PRKMK1