



ELK Biotechnology

p120 Rabbit pAb
Catalog Number: ES8534
For research use only.

Overview

Product Name	p120 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthetic peptide from human protein at AA range; 841-890
Specificity	The antibody detects endogenous p120 protein
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	p120
Gene Name	CTNND1 KIAA0384
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	108kDa
Human Gene ID	1500
Human Swiss-Prot Number	O60716
Alternative Names	CTNND1 KIAA0384