



ELK Biotechnology

HSC70 Rabbit Rabbit pAb
Catalog Number: ES8558
For research use only.

Overview

Product Name	HSC70 Rabbit polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/50-1/300
Immunogen	Recombinant Protein of Heat shock cognate 71 kDa protein (Heat shock 70 kDa protein 8)
Specificity	The antibody detects endogenous HSC70 protein
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	HSC70 Rabbit
Gene Name	HSPA8 HSC70 HSP73 HSPA10
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	70-72kDa
Human Gene ID	3312
Human Swiss-Prot Number	P11142
Alternative Names	Heat shock cognate 71 kDa protein (Heat shock 70 kDa protein 8)