



ELK Biotechnology

JNK2 Rabbit Rabbit pAb
Catalog Number: ES8564
For research use only.

Overview

Product Name	JNK2 Rabbit polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/50-1/300
Immunogen	Recombinant Protein of JNK2
Specificity	The antibody detects endogenous JNK2 protein
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	JNK2 Rabbit
Gene Name	MAPK9 JNK2 PRKM9 SAPK1A
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	46,54kDa
Human Gene ID	5601
Human Swiss-Prot Number	P45984
Alternative Names	Mitogen-activated protein kinase 9 (MAP kinase 9) (MAPK 9) (EC 2.7.11.24) (JNK-55) (Stress-activated protein kinase 1a) (SAPK1a) (Stress-activated protein kinase JNK2) (c-Jun N-terminal kinase 2)