

Nrf2 Rabbit Rabbit pAb Catalog Number: ES8568 For research use only.

Overview

Product Name	Nrf2 Rabbit polyclonal antibody
Host species	Rabbit
Applications	ІНС-р
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/50-1/300
Immunogen	Recombinant Protein of Nrf2
Specificity	The antibody detects endogenous Nrf2 protein
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Nrf2 Rabbit
Gene Name	NFE2L2 NRF2
Purification	The antibody was affinity-purified from rabbit serum by affinity- chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	75-100kDa
Human Gene ID	4780
Human Swiss-Prot Number	Q16236
Alternative Names	Nuclear factor erythroid 2-related factor 2 (NF-E2-related factor 2) (NFE2-related factor 2) (HEBP1) (Nuclear factor, erythroid derived 2, like 2)