



ELK Biotechnology

ATG4a Rabbit Rabbit pAb
Catalog Number: ES8569
For research use only.

Overview

Product Name	ATG4a Rabbit polyclonal antibody
Host species	Rabbit
Applications	IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/50-1/300
Immunogen	Recombinant Protein of ATG4a
Specificity	The antibody detects endogenous ATG4a protein
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ATG4a Rabbit
Gene Name	ATG4A APG4A AUTL2
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	45kDa
Human Gene ID	115201
Human Swiss-Prot Number	Q8WYN0
Alternative Names	Cysteine protease ATG4A (EC 3.4.22.-) (AUT-like 2 cysteine endopeptidase) (Autophagin-2) (Autophagy-related cysteine endopeptidase 2) (Autophagy-related protein 4 homolog A) (hAPG4A)