



ELK Biotechnology

ATG4b Rabbit Rabbit pAb
Catalog Number: ES8570
For research use only.

Overview

Product Name	ATG4b Rabbit polyclonal antibody
Host species	Rabbit
Applications	IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/50-1/300
Immunogen	Recombinant Protein of ATG4b
Specificity	The antibody detects endogenous ATG4b protein
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ATG4b Rabbit
Gene Name	ATG4B APG4B AUTL1 KIAA0943
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	44kDa
Human Gene ID	23192
Human Swiss-Prot Number	Q9Y4P1
Alternative Names	Cysteine protease ATG4B (EC 3.4.22.-) (AUT-like 1 cysteine endopeptidase) (Autophagin-1) (Autophagy-related cysteine endopeptidase 1) (Autophagy-related protein 4 homolog B) (hAPG4B)