



ELK Biotechnology

ATG7 Rabbit Rabbit pAb
Catalog Number: ES8573
For research use only.

Overview

Product Name	ATG7 Rabbit polyclonal antibody
Host species	Rabbit
Applications	IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/50-1/300
Immunogen	Recombinant Protein of ATG7
Specificity	The antibody detects endogenous ATG7 protein
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ATG7 Rabbit
Gene Name	ATG7 APG7L
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	78kDa
Human Gene ID	10533
Human Swiss-Prot Number	O95352
Alternative Names	Ubiquitin-like modifier-activating enzyme ATG7 (ATG12-activating enzyme E1 ATG7) (Autophagy-related protein 7) (APG7-like) (hAGP7) (Ubiquitin-activating enzyme E1-like protein)