



# ELK Biotechnology

**ATG14L Rabbit Rabbit pAb**  
**Catalog Number: ES8576**  
**For research use only.**

## Overview

---

<b>Product Name</b>	ATG14L Rabbit polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>IHC-p</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Immunohistochemistry:1/50-1/300</b>
<b>Immunogen</b>	Recombinant Protein of ATG14L
<b>Specificity</b>	The antibody detects endogenous ATG14L protein
<b>Formulation</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	ATG14L Rabbit
<b>Gene Name</b>	ATG14 KIAA0831
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>55kDa</b>
<b>Human Gene ID</b>	22863
<b>Human Swiss-Prot Number</b>	Q6ZNE5
<b>Alternative Names</b>	Beclin 1-associated autophagy-related key regulator (Barkor) (Autophagy-related protein 14-like protein) (Atg14L)