



# ELK Biotechnology

**c-Abl Rabbit Rabbit pAb**  
**Catalog Number: ES8577**  
**For research use only.**

## Overview

---

<b>Product Name</b>	c-Abl Rabbit polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>IHC-p</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Immunohistochemistry:1/50-1/300</b>
<b>Immunogen</b>	Recombinant Protein of c-Abl
<b>Specificity</b>	The antibody detects endogenous c-Abl protein
<b>Formulation</b>	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	c-Abl Rabbit
<b>Gene Name</b>	ABL1 ABL JTK7
<b>Purification</b>	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>120kDa</b>
<b>Human Gene ID</b>	25
<b>Human Swiss-Prot Number</b>	P00519
<b>Alternative Names</b>	Tyrosine-protein kinase ABL1 (EC 2.7.10.2) (Abelson murine leukemia viral oncogene homolog 1) (Abelson tyrosine-protein kinase 1) (Proto-oncogene c-Abl) (p150)