



ELK Biotechnology

NBR1 Rabbit Rabbit pAb
Catalog Number: ES8578
For research use only.

Overview

Product Name	NBR1 Rabbit polyclonal antibody
Host species	Rabbit
Applications	IHC-p
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/50-1/300
Immunogen	Recombinant Protein of NBR1
Specificity	The antibody detects endogenous NBR1 protein
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	NBR1 Rabbit
Gene Name	NBR1 1A13B KIAA0049 M17S2 MIG19
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	120kDa
Human Gene ID	4077
Human Swiss-Prot Number	Q14596
Alternative Names	Next to BRCA1 gene 1 protein (Cell migration-inducing gene 19 protein) (Membrane component chromosome 17 surface marker 2) (Neighbor of BRCA1 gene 1 protein) (Protein 1A1-3B)