



ELK Biotechnology

IRF3 Rabbit pAb
Catalog Number: ES8598
For research use only.

Overview

Product Name	IRF3 polyclonal antibody
Host species	Rabbit
Applications	WB,IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 Immunohistochemistry:1/500-1/200 ELISA 1/10000-1/20000
Immunogen	Synthetic peptide from human protein at AA range; 1268-1341
Specificity	The antibody detects endogenous IRF3
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	IRF3
Gene Name	IRF3
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	3661
Human Swiss-Prot Number	Q14653
Alternative Names	Interferon regulatory factor 3 (IRF-3)