



## ELK Biotechnology

**Catenin- $\beta$  Rabbit pAb**  
**Catalog Number: ES8621**  
**For research use only.**

### Overview

---

Product Name	Catenin- $\beta$ polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthetic peptide from human protein at AA range; 20-80
Specificity	The antibody detects endogenous Catenin- $\beta$
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Catenin- $\beta$
Gene Name	CTNNB1 CTNNB OK/SW-cl.35 PRO2286
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	90kDa
Human Gene ID	1499
Human Swiss-Prot Number	P35222
Alternative Names	Catenin beta-1 (Beta-catenin)