



ELK Biotechnology

Histone H2B Rabbit pAb
Catalog Number: ES8623
For research use only.

Overview

Product Name	Histone H2B polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthetic peptide from human protein at AA range; 1-50
Specificity	The antibody detects endogenous Histone H2B
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Histone H2B
Gene Name	Histone H2B
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	15kDa
Human Gene ID	255626
Human Swiss-Prot Number	Q96A08/P33778/P62807
Alternative Names	