



ELK Biotechnology

IFN- α / β R α Rabbit pAb
Catalog Number: ES8642
For research use only.

Overview

Product Name	IFN- α / β R α polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Immunohistochemistry:1/50-1/200 ELISA:1/10000-1/20000
Immunogen	Synthetic peptide from human protein at AA range; 300-350
Specificity	The antibody detects endogenous IFN- α / β R α
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	IFN- α / β R α
Gene Name	IFNAR1 IFNAR
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	3454
Human Swiss-Prot Number	P17181
Alternative Names	Interferon alpha/beta receptor 1 (IFN-R-1) (IFN-alpha/beta receptor 1) (Cytokine receptor class-II member 1) (Cytokine receptor family 2 member 1) (CRF2-1) (Type I interferon receptor 1)