



ELK Biotechnology

GRO α Rabbit pAb
Catalog Number: ES8668
For research use only.

Overview

Product Name	GRO α polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Immunohistochemistry:1/50-1/200 ELISA:1/10000-1/20000
Immunogen	Synthetic peptide from human protein at AA range; 70-107
Specificity	The antibody detects endogenous GRO α
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GRO α
Gene Name	CXCL1 GRO GRO1 GROA MGSA SCYB1
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	2919
Human Swiss-Prot Number	P09341
Alternative Names	Growth-regulated alpha protein (C-X-C motif chemokine 1) (GRO-alpha(1-73)) (Melanoma growth stimulatory activity) (MGSA) (Neutrophil-activating protein 3) (NAP-3) [Cleaved into; GRO-alpha(4-73); GRO-alpha(5-73); GRO-alpha(6-73)]