



ELK Biotechnology

ACE1 Rabbit pAb
Catalog Number: ES8680
For research use only.

Overview

Product Name	ACE1 polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse
Recommended dilutions	Immunohistochemistry:1/50-1/200 ELISA:1/10000-1/20000
Immunogen	Synthetic peptide from human protein at AA range; 41-90
Specificity	The antibody detects endogenous ACE1
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ACE1
Gene Name	ACE DCP DCP1
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	1636
Human Swiss-Prot Number	P12821
Alternative Names	Angiotensin-converting enzyme (ACE) (EC 3.2.1.-) (EC 3.4.15.1) (Dipeptidyl carboxypeptidase I) (Kininase II) (CD antigen CD143) [Cleaved into; Angiotensin-converting enzyme, soluble form]