



# ELK Biotechnology

**ZAP-70 Rabbit pAb**  
**Catalog Number: ES8738**  
**For research use only.**

## Overview

---

Product Name	ZAP-70 polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	<b>Immunohistochemistry:1/50-1/200</b> <b>ELISA:1/10000-1/20000</b>
Immunogen	Synthetic peptide from human protein at AA range; 391-440
Specificity	The antibody detects endogenous ZAP-70
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	ZAP-70
Gene Name	ZAP70 SRK
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	<b>kDa</b>
Human Gene ID	7535
Human Swiss-Prot Number	P43403
Alternative Names	Tyrosine-protein kinase ZAP-70 (EC 2.7.10.2) (70 kDa zeta-chain associated protein) (Syk-related tyrosine kinase)