



ELK Biotechnology

Flt-4 Rabbit pAb
Catalog Number: ES8770
For research use only.

Overview

Product Name	Flt-4 polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/50-1/200 ELISA:1/10000-1/20000
Immunogen	Synthetic peptide from human protein at AA range; 720-760
Specificity	The antibody detects endogenous Flt-4
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Flt-4
Gene Name	FLT4 VEGFR3
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	2324
Human Swiss-Prot Number	P35916
Alternative Names	Vascular endothelial growth factor receptor 3 (VEGFR-3) (EC 2.7.10.1) (Fms-like tyrosine kinase 4) (FLT-4) (Tyrosine-protein kinase receptor FLT4)