



ELK Biotechnology

CUL-1 Rabbit pAb
Catalog Number: ES8776
For research use only.

Overview

Product Name	CUL-1 polyclonal antibody
Host species	Rabbit
Applications	IHC-p,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Immunohistochemistry:1/50-1/200 ELISA:1/10000-1/20000
Immunogen	Synthetic peptide from human protein at AA range; 130-180
Specificity	The antibody detects endogenous CUL-1
Formulation	PBS, pH 7.4, containing 0.02% sodium azide as Preservative and 50% Glycerol.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	CUL-1
Gene Name	CUL1
Purification	The antibody was affinity-purified from rabbit serum by affinity-chromatography using specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	8454
Human Swiss-Prot Number	Q13616
Alternative Names	Cullin-1 (CUL-1)