



ELK Biotechnology

CD22 Rabbit pAb
Catalog Number: ES8813
For research use only.

Overview

Product Name	CD22 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/1000
Immunogen	Synthesized peptide derived from human CD22. at AA range; 780-860
Specificity	This antibody detects endogenous levels of CD22
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	CD22
Gene Name	CD22 SIGLEC2
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	933
Human Swiss-Prot Number	P20273
Alternative Names	B-cell receptor CD22 (B-lymphocyte cell adhesion molecule) (BL-CAM) (Sialic acid-binding Ig-like lectin 2) (Siglec-2) (T-cell surface antigen Leu-14) (CD antigen CD22)