



# ELK Biotechnology

**E2F-1 Rabbit pAb**  
**Catalog Number: ES8821**  
**For research use only.**

## Overview

---

<b>Product Name</b>	E2F-1 polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB,ELISA</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/1000</b>
<b>Immunogen</b>	Synthesized peptide derived from human E2F-1. at AA range; 100-140
<b>Specificity</b>	This antibody detects endogenous levels of E2F-1
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	E2F-1
<b>Gene Name</b>	E2F1 RBBP3
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>kDa</b>
<b>Human Gene ID</b>	1869
<b>Human Swiss-Prot Number</b>	Q01094
<b>Alternative Names</b>	Transcription factor E2F1 (E2F-1) (PBR3) (Retinoblastoma-associated protein 1) (RBAP-1) (Retinoblastoma-binding protein 3) (RBBP-3) (pRB-binding protein E2F-1)