



ELK Biotechnology

FoxO3A Rabbit pAb
Catalog Number: ES8823
For research use only.

Overview

Product Name	FoxO3A polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/1000 icc:1/50
Immunogen	Synthesized peptide derived from human FoxO3A. at AA range; 550-580
Specificity	This antibody detects endogenous levels of FoxO3A
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	FoxO3A
Gene Name	FOXO3 FKHL1 FOXO3A
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	2309
Human Swiss-Prot Number	O43524
Alternative Names	Forkhead box protein O3 (AF6q21 protein) (Forkhead in rhabdomyosarcoma-like 1)