



ELK Biotechnology

FoxO4 Rabbit pAb
Catalog Number: ES8824
For research use only.

Overview

Product Name	FoxO4 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/1000 icc:1/400
Immunogen	Synthesized peptide derived from human FoxO4. at AA range; 170-200
Specificity	This antibody detects endogenous levels of FoxO4
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	FoxO4
Gene Name	FOXO4 AFX AFX1 MLLT7
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	4303
Human Swiss-Prot Number	P98177
Alternative Names	Forkhead box protein O4 (Fork head domain transcription factor AFX1)