



ELK Biotechnology

GI Syn Rabbit pAb
Catalog Number: ES8835
For research use only.

Overview

Product Name	GI Syn polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/50000
Immunogen	Synthesized peptide derived from human GI Syn. at AA range; 301-350
Specificity	This antibody detects endogenous levels of GI Syn
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	GI Syn
Gene Name	GLUL GLNS
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	kDa
Human Gene ID	2752
Human Swiss-Prot Number	P15104
Alternative Names	Glutamine synthetase (GS) (EC 6.3.1.2) (Glutamate decarboxylase) (EC 4.1.1.15) (Glutamate--ammonia ligase)