



ELK Biotechnology

MMP-10 Rabbit pAb
Catalog Number: ES8852
For research use only.

Overview

Product Name	MMP-10 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthesized peptide derived from human MMP-10. at AA range; 351-400
Specificity	MMP-10 Polyclonal Antibody detects endogenous levels of MMP-10
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	MMP-10
Gene Name	MMP10
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	52kDa
Human Gene ID	4319
Human Swiss-Prot Number	P09238
Alternative Names	Stromelysin-2 (SL-2) (EC 3.4.24.22) (Matrix metalloproteinase-10) (MMP-10) (Transin-2)