



ELK Biotechnology

Cdk5 Rabbit pAb
Catalog Number: ES8853
For research use only.

Overview

Product Name	Cdk5 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthesized peptide derived from human Cdk5. at AA range; 1-50
Specificity	Cdk5 Polyclonal Antibody detects endogenous levels of Cdk5
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Cdk5
Gene Name	CDK5
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	33kDa
Human Gene ID	1020
Human Swiss-Prot Number	Q00535
Alternative Names	Cyclin-dependent kinase 5 (EC 2.7.11.22) (Cell division protein kinase 5) (Serine/threonine-protein kinase PSSALRE) (Tau protein kinase II catalytic subunit) (TPKII catalytic subunit)