



ELK Biotechnology

COL19A1 Rabbit pAb
Catalog Number: ES8856
For research use only.

Overview

Product Name	COL19A1 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthesized peptide derived from human COL19A1. at AA range; 161-210
Specificity	COL19A1 Polyclonal Antibody detects endogenous levels of COL19A1
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	COL19A1
Gene Name	COL19A1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	130kDa
Human Gene ID	1310
Human Swiss-Prot Number	Q14993
Alternative Names	Collagen alpha-1(XIX) chain (Collagen alpha-1(Y) chain)