



ELK Biotechnology

Maspin Rabbit pAb
Catalog Number: ES8868
For research use only.

Overview

Product Name	Maspin polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthesized peptide derived from human Maspin. at AA range; 201-250
Specificity	Maspin Polyclonal Antibody detects endogenous levels of Maspin
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Maspin
Gene Name	SERPINB5
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	42kDa
Human Gene ID	5268
Human Swiss-Prot Number	P36952
Alternative Names	Serpin B5 (Maspin) (Peptidase inhibitor 5) (PI-5)