



ELK Biotechnology

p53 Rabbit pAb
Catalog Number: ES8902
For research use only.

Overview

| | |
|---------------------------------|---|
| Product Name | p53 polyclonal antibody |
| Host species | Rabbit |
| Applications | WB,ELISA |
| Species Cross-Reactivity | Human,Mouse,Rat |
| Recommended dilutions | Western Blot:1/500-1/2000 ELISA 1/10000-1/20000 |
| Immunogen | Synthesized peptide derived from human p53 Polyclonal |
| Specificity | This antibody detects endogenous levels of p53. |
| Formulation | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide. |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| Protein Name | p53 |
| Gene Name | TP53 P53 |
| Purification | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| Clonality | Polyclonal |
| Concentration | 1 mg/ml |
| Observed band | 53kDa |
| Human Gene ID | 7157 |
| Human Swiss-Prot Number | P04637 |
| Alternative Names | Cellular tumor antigen p53 (Antigen NY-CO-13) (Phosphoprotein p53) (Tumor suppressor p53) |