



# ELK Biotechnology

**Src Rabbit pAb**  
**Catalog Number: ES8903**  
**For research use only.**

## Overview

---

<b>Product Name</b>	Src polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB,ELISA</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/500-1/2000</b> <b>ELISA</b> <b>1/10000-1/20000</b>
<b>Immunogen</b>	Synthesized peptide derived from human Src Polyclonal
<b>Specificity</b>	This antibody detects endogenous levels of Src.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	Src
<b>Gene Name</b>	SRC SRC1
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>60kDa</b>
<b>Human Gene ID</b>	6714
<b>Human Swiss-Prot Number</b>	P12931
<b>Alternative Names</b>	Proto-oncogene tyrosine-protein kinase Src (EC 2.7.10.2) (Proto-oncogene c-Src) (pp60c-src) (p60-Src)