



ELK Biotechnology

Bcl-6 Rabbit pAb
Catalog Number: ES8934
For research use only.

Overview

Product Name	Bcl-6 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthesized peptide derived from B-cell lymphoma 6 protein at AA range; 261-310
Specificity	This antibody detects endogenous levels of Bcl-6.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	Bcl-6
Gene Name	BCL6
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	90kDa
Human Gene ID	604
Human Swiss-Prot Number	P41182
Alternative Names	