



ELK Biotechnology

PAR4 Rabbit pAb
Catalog Number: ES8940
For research use only.

Overview

Product Name	PAR4 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/10000-1/20000
Immunogen	Synthesized peptide derived from human PAR4 Polyclonal
Specificity	This antibody detects endogenous levels of PAR4.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	PAR4
Gene Name	F2RL3 PAR4
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	40kDa
Human Gene ID	9002
Human Swiss-Prot Number	Q96RI0
Alternative Names	Proteinase-activated receptor 4 (PAR-4) (Coagulation factor II receptor-like 3) (Thrombin receptor-like 3)