



# ELK Biotechnology

**Caspase 3 p17 Rabbit pAb**  
**Catalog Number: ES8941**  
**For research use only.**

## Overview

---

<b>Product Name</b>	Caspase 3 p17 polyclonal antibody
<b>Host species</b>	Rabbit
<b>Applications</b>	<b>WB,ELISA</b>
<b>Species Cross-Reactivity</b>	<b>Human,Mouse,Rat</b>
<b>Recommended dilutions</b>	<b>Western Blot:1/500-1/2000</b> <b>ELISA</b> <b>1/10000-1/20000</b>
<b>Immunogen</b>	Synthesized peptide derived from human Caspase 3 p17 Polyclonal
<b>Specificity</b>	This antibody detects endogenous levels of Caspase 3 p17.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
<b>Protein Name</b>	Caspase 3 p17
<b>Gene Name</b>	CASP3 CPP32
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	<b>17kDa</b>
<b>Human Gene ID</b>	836
<b>Human Swiss-Prot Number</b>	P42574
<b>Alternative Names</b>	Caspase-3 (CASP-3) (EC 3.4.22.56) (Apopain) (Cysteine protease CPP32) (CPP-32) (Protein Yama) (SREBP cleavage activity 1) (SCA-1) [Cleaved into; Caspase-3 subunit p17; Caspase-3 subunit p12]