



ELK Biotechnology

FGFR1 Rabbit pAb
Catalog Number: ES8963
For research use only.

Overview

Product Name	FGFR1 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human,Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/5000-1/20000
Immunogen	Synthesized peptide derived from human protein . at AA range; 600-680
Specificity	FGFR1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	FGFR1
Gene Name	FGFR1 BFGFR CEK FGFBR FLG FLT2 HBGFR
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	90kDa
Human Gene ID	2260
Human Swiss-Prot Number	P11362
Alternative Names	