



ELK Biotechnology

TRPV1 Rabbit pAb
Catalog Number: ES9015
For research use only.

Overview

Product Name	TRPV1 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Mouse,Rat
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/5000-1/20000
Immunogen	Synthesized peptide derived from human protein . at AA range; 440-520
Specificity	TRPV1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	TRPV1
Gene Name	TRPV1 VR1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	92kDa
Human Gene ID	7442
Human Swiss-Prot Number	Q8NER1
Alternative Names	