



ELK Biotechnology

WIT1 Rabbit pAb

Catalog Number: ES9049

For research use only.

Overview

Product Name	WIT1 polyclonal antibody
Host species	Rabbit
Applications	WB,ELISA
Species Cross-Reactivity	Human
Recommended dilutions	Western Blot:1/500-1/2000 ELISA 1/5000-1/20000
Immunogen	Synthesized peptide derived from human protein . at AA range; 20-100
Specificity	WIT1 Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, and 0.02% sodium azide.
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
Protein Name	WIT1
Gene Name	WT1-AS WIT1
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	10kDa
Human Gene ID	0
Human Swiss-Prot Number	Q06250
Alternative Names	