



# FAKD5 rabbit pAb

Cat No.:ES9075

For research use only

## Overview

<b>Product Name</b>	FAKD5 rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB;ELISA
<b>Species Cross-Reactivity</b>	Human;Mouse
<b>Recommended dilutions</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Immunogen</b>	Synthesized peptide derived from human protein . at AA range: 500-580
<b>Specificity</b>	FAKD5 Polyclonal Antibody detects endogenous levels of protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C . Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	FAST kinase domain-containing protein 5
<b>Gene Name</b>	FASTKD5 KIAA1792
<b>Cellular localization</b>	Mitochondrion matrix, mitochondrion nucleoid . Localizes to mitochondrial RNA granules found in close proximity to the mitochondrial nucleoids. .
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	84kD
<b>Human Gene ID</b>	60493
<b>Human Swiss-Prot Number</b>	Q7L8L6
<b>Alternative Names</b>	
<b>Background</b>	similarity:Belongs to the FAST kinase family.,similarity:Contains 1 RAP domain.,

