

IFN17 rabbit pAb

Cat No.: ES9164

For research use only

Overview

Product Name IFN17 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human;Rat;Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein .

at AA range: 1-80

Specificity IFN17 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

StorageStore at -20°C. Avoid repeated freeze-thaw cycles.Protein NameInterferon alpha-17 (IFN-alpha-17) (Interferon

alpha-88) (Interferon alpha-I') (LeIF I) (Interferon

alpha-T)

Gene Name IFNA17 **Cellular localization** Secreted.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 20kD
Human Gene ID 3451
Human Swiss-Prot Number P01571

Alternative Names

Background function: Produced by macrophages, IFN-alpha have

antiviral activities. Interferon stimulates the

production of two enzymes: a protein kinase and an oligoadenylate synthetase., similarity: Belongs to the


alpha/beta interferon family.,



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com





Western blot analysis of lysates from SH-SY5Y cells, primary antibody was diluted at 1:1000, 4° over night

