

PDC10 rabbit pAb

Cat No.:ES9217

For research use only

Overview

Product Name PDC10 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein .

at AA range: 40-120

Specificity PDC10 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

StorageStore at -20°C. Avoid repeated freeze-thaw cycles.Protein NameProgrammed cell death protein 10 (Cerebral

cavernous malformations 3 protein) (TF-1 cell

apoptosis-related protein 15)

Gene Name PDCD10 CCM3 TFAR15

Cellular localization Cytoplasm. Golgi apparatus membrane; Peripheral

membrane protein; Cytoplasmic side. Cell membrane; Peripheral membrane protein; Cytoplasmic side. Partially co-localizes with

endogenous PXN at the leading edges of migrating

cells.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/mlObserved band23kDHuman Gene ID11235Human Swiss-Prot NumberQ9BUL8

Alternative Names

Background This gene encodes an evolutionarily conserved


protein associated with cell apoptosis. The protein





interacts with the serine/threonine protein kinase MST4 to modulate the extracellular signal-regulated kinase (ERK) pathway. It also interacts with and is phosphoryated by serine/threonine kinase 25, and is thought to function in a signaling pathway essential for vascular developent. Mutations in this gene are one cause of cerebral cavernous malformations, which are vascular malformations that cause seizures and cerebral hemorrhages. Multiple alternatively spliced variants, encoding the same protein, have been identified. [provided by RefSeq, Jul 2008],

Western blot analysis of lysates from SW480 cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950

