

SYT2 rabbit pAb

Cat No.: ES9241

For research use only

Overview

Product Name SYT2 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein .

at AA range: 40-120

Specificity SYT2 Polyclonal Antibody detects endogenous levels

of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Synaptotagmin-2 (Synaptotagmin II) (SytII)

Gene Name SYT2

Cellular localization Cytoplasmic vesicle, secretory vesicle, synaptic

vesicle membrane ; Single-pass membrane protein . Cytoplasmic vesicle, secretory vesicle, chromaffin granule membrane ; Single-pass membrane protein .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

ClonalityPolyclonalConcentration1 mg/mlObserved band46kDHuman Gene ID127833Human Swiss-Prot NumberQ8N9I0

Alternative Names

Background This gene encodes a synaptic vesicle membrane

protein. The encoded protein is thought to function as a calcium sensor in vesicular trafficking and exocytosis. Mutations in this gene are associated with myasthenic syndrome, presynaptic, congenital,

with or without motor neuropathy. Alternative



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com



MDa 1
180 -140 -100 -75 -60 -45 -25 -15 -10 --

splicing results in multiple transcript variants. [provided by RefSeq, Oct 2014],

Western blot analysis of lysates from Jarkat cells, primary antibody was diluted at 1:1000, 4° over night



+86-27-59760950