

RM02 rabbit pAb

Cat No.: ES9299

For research use only

Overview

Product Name RM02 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein . at

AA range: 210-290

Specificity RM02 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name 39S ribosomal protein L2, mitochondrial (L2mt)

(MRP-L2)

Gene Name MRPL2 CGI-22
Cellular localization Mitochondrion .

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 33kD
Human Gene ID 51069
Human Swiss-Prot Number Q5T653

Alternative Names

Background Mammalian mitochondrial ribosomal proteins are

encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition

compared to prokaryotic ribosomes, where this ratio is reversed. Another difference between mammalian



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com



mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein that belongs to the EcoL2 ribosomal protein family. A pseudogene corresponding to this gene is found on chromosome 12q. Alternative splicing results in multiple transcript variants encoding distinct isoforms.



+86-27-59760950