

RM43 rabbit pAb

Cat No.:ES9304

For research use only

Overview

Product Name RM43 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Rat; Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein.

at AA range: 120-200

Specificity RM43 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

StorageStore at -20°C. Avoid repeated freeze-thaw cycles.Protein Name39S ribosomal protein L43, mitochondrial (L43mt)

(MRP-L43) (Mitochondrial ribosomal protein

bMRP36a)

Gene Name MRPL43

Cellular localization Mitochondrion.

Purification The antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 23kD
Human Gene ID 84545
Human Swiss-Prot Number Q8N983

Alternative Names

Background Mammalian mitochondrial ribosomal proteins are

encoded by nuclear genes and help in protein synthesis within the mitochondrion. Mitochondrial ribosomes (mitoribosomes) consist of a small 28S subunit and a large 39S subunit. They have an estimated 75% protein to rRNA composition

compared to prokaryotic ribosomes, where this ratio



+86-27-59760950


ELKbio@ELKbiotech.com

www.elkbiotech.com



is reversed. Another difference between mammalian mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. This gene and the gene for a semaphorin class 4 protein (SEMA4G) overlap at map location 10q24.31 and are transcribed in opposite directions. Sequence analysis identified multiple transcript variants encodin

Western blot analysis of lysates from PC12 cells, primary antibody was diluted at 1:1000, 4° over night





+86-27-59760950