

RM55 rabbit pAb

Cat No.: ES9307

For research use only

Overview

Product Name	RM55 rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Mouse
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from human protein .
-	at AA range: 30-110
Specificity	RM55 Polyclonal Antibody detects endogenous
	levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and
	0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	39S ribosomal protein L55, mitochondrial (L55mt)
	(MRP-L55)
Gene Name	MRPL55 UNQ5835/PRO19675
Cellular localization	Mitochondrion .
Purification	The antibody was affinity-purified from rabbit
	antiserum by affinity-chromatography using
	epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	14kD
Human Gene ID	128308
Human Swiss-Prot Number	Q7Z7F7
Alternative Names	
Background	Mammalian mitochondrial ribosomal proteins are
	encoded by nuclear genes and help in protein
	synthesis within the mitochondrion. Mitochondrial
	ribosomes (mitoribosomes) consist of a small 28S
	subunit and a large 39S subunit. They have an
	estimated 75% protein to rRNA composition
	compared to prokaryotic ribosomes, where this ratio
	· · · · · · · · · · · · · · · · · · ·



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com


is reversed. Another difference between mammalian

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



mitoribosomes and prokaryotic ribosomes is that the latter contain a 5S rRNA. Among different species, the proteins comprising the mitoribosome differ greatly in sequence, and sometimes in biochemical properties, which prevents easy recognition by sequence homology. This gene encodes a 39S subunit protein. Multiple transcript variants encoding two different isoforms were identified through sequence analysis. [provided by RefSeq, Jul 2008],

Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night





+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C