

ABHD1 rabbit pAb

Cat No.: ES9333

For research use only

Overview

Product Name ABHD1 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human; Mouse; Rat

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from human protein .

at AA range: 220-300

Specificity ABHD1 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

Storage Store at -20°C. Avoid repeated freeze-thaw cycles.

Protein Name Abhydrolase domain-containing protein 1 (EC

3.1.1.-) (Lung alpha/beta hydrolase 1)

Gene Name ABHD1 LABH1

Cellular localizationMembrane ; Single-pass type II membrane protein .PurificationThe antibody was affinity-purified from rabbit antiserum by affinity-chromatography using

and serum by annity-chromatography

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 44kD
Human Gene ID 84696
Human Swiss-Prot Number Q96SE0

Alternative Names

Background This gene is a member of the AB hydrolase

superfamily and encodes a protein with an


alpha/beta hydrolase fold. This domain is common to a number of hydrolytic enzymes of widely differing phylogenetic origins and catalytic functions. [provided by RefSeq, Jul 2008],



+86-27-59760950 ELKbio@ELKbiotech.com

www.elkbiotech.com





Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

