

CEA19 rabbit pAb

Cat No.: ES9505

For research use only

Overview

Product Name CEA19 rabbit pAb

Host species Rabbit
Applications WB;ELISA

Species Cross-Reactivity Human;Rat;Mouse;

Recommended dilutions WB 1:500-2000 ELISA 1:5000-20000

Immunogen Synthesized peptide derived from part region of

human protein

Specificity CEA19 Polyclonal Antibody detects endogenous

levels of protein.

Formulation Liquid in PBS containing 50% glycerol, 0.5% BSA and

0.02% sodium azide.

StorageStore at -20°C. Avoid repeated freeze-thaw cycles.Protein NameCarcinoembryonic antigen-related cell adhesion

molecule 19 (Carcinoembryonic antigen-like 1)

Gene Name CEACAM19 CEAL1 UNQ2973/PRO7436

Cellular localizationMembrane ; Single-pass type I membrane protein .PurificationThe antibody was affinity-purified from rabbit

antiserum by affinity-chromatography using

epitope-specific immunogen.

Clonality Polyclonal
Concentration 1 mg/ml
Observed band 33kD
Human Gene ID 56971
Human Swiss-Prot Number Q7Z692

Alternative Names

Background induction: Down-regulated by dexamethasone (in

vitro).,similarity:Belongs to the immunoglobulin superfamily. CEA family.,tissue specificity:Ubiquitous with highest expression in prostate, uterus, fetal brain, mammary gland, adrenal gland, skeletal muscle, small intestine, and kidney, and lower expression in lung, cerebellum, testis, liver,

pancreas, bone marrow and ovary.,




+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com





Western blot analysis of lysates from K562 cells, primary antibody was diluted at 1:1000, 4° over night

