

CTSR4 rabbit pAb

Cat No.:ES9512

For research use only

Overview

Product Name	CTSR4 rabbit pAb	
Host species	Rabbit	
Applications	WB;ELISA	
Species Cross-Reactivity	Human;Rat;Mouse;	
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000	
Immunogen	Synthesized peptide derived from part region of	
-	human protein	
Specificity	CTSR4 Polyclonal Antibody detects endogenous	
	levels of protein.	
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and	
	0.02% sodium azide.	
Storage	Store at -20 $^\circ C$. Avoid repeated freeze-thaw cycles.	
Protein Name	Cation channel sperm-associated protein 4	
	(CatSper4)	
Gene Name	CATSPER4	
Cellular localization	Cell projection, cilium, flagellum membrane ;	
	Multi-pass membrane protein . Specifically located	
	in the principal piece of the sperm tail	
Purification	The antibody was affinity-purified from rabbit	
	antiserum by affinity-chromatography using	
	epitope-specific immunogen.	
Clonality	Polyclonal	
Concentration	1 mg/ml	
Observed band	51kD	
Human Gene ID	378807	
Human Swiss-Prot Number	Q7RTX7	
Alternative Names		
Background	function:Voltage-gated calcium channel that plays a	
	central role in sperm cell hyperactivation. Controls	
	calcium entry to mediate the hyperactivated	
	motility, a step needed for sperm motility which is	
	essential late in the preparation of sperm for	
	fertilization. Activated by intracellular	



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C



alkalinization.,similarity:Belongs to the cation channel sperm-associated family.,subunit:Heterotetramer; possibly composed of CATSPER1, CATSPER2, CATSPER3 and CATSPER4 (Potential). Interacts with CATSPER1.,tissue specificity:Testis-specific. Specifically present in the principal piece of sperm tail.,



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road, Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C