



GON4L rabbit pAb

Cat No.:ES9679

For research use only

Overview

Product Name	GON4L rabbit pAb
Host species	Rabbit
Applications	WB;ELISA
Species Cross-Reactivity	Human;Rat;Mouse;
Recommended dilutions	WB 1:500-2000 ELISA 1:5000-20000
Immunogen	Synthesized peptide derived from part region of human protein
Specificity	GON4L Polyclonal Antibody detects endogenous levels of protein.
Formulation	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Storage	Store at -20°C. Avoid repeated freeze-thaw cycles.
Protein Name	GON-4-like protein (GON-4 homolog)
Gene Name	GON4L GON4 KIAA1606
Cellular localization	Nucleus .
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Clonality	Polyclonal
Concentration	1 mg/ml
Observed band	246kD
Human Gene ID	54856
Human Swiss-Prot Number	Q3T8J9
Alternative Names	
Background	function:Enhances transcription activation by YY1.,sequence caution:Translated as Tyr.,similarity:Contains 1 Myb-like domain.,similarity:Contains 2 PAH (paired amphipathic helix) repeats.,subunit:Interacts with YY1.,tissue specificity:Ubiquitous. Detected in small intestine, skeletal muscle, lung, pancreas, brain, stomach, spleen, colon and heart. Detected at very low levels in healthy liver. Highly expressed in most





ELK Biotechnology

liver carcinomas.,



+86-27-59760950

ELKbio@ELKbiotech.com

www.elkbiotech.com

23-2, No.388 Gaoxin 2nd Road,Wuhan East Lake Hi-tech Development Zone, Hubei , P.R.C