



# NOC4L rabbit pAb

Cat No.:ES9950

For research use only

## Overview

<b>Product Name</b>	NOC4L rabbit pAb
<b>Host species</b>	Rabbit
<b>Applications</b>	WB;ELISA
<b>Species Cross-Reactivity</b>	Human;Rat;Mouse;
<b>Recommended dilutions</b>	WB 1:500-2000 ELISA 1:5000-20000
<b>Immunogen</b>	Synthesized peptide derived from human protein . at AA range: 370-450
<b>Specificity</b>	NOC4L Polyclonal Antibody detects endogenous levels of protein.
<b>Formulation</b>	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
<b>Storage</b>	Store at -20°C . Avoid repeated freeze-thaw cycles.
<b>Protein Name</b>	Nucleolar complex protein 4 homolog (NOC4 protein homolog) (NOC4-like protein) (Nucleolar complex-associated protein 4-like protein)
<b>Gene Name</b>	NOC4L
<b>Cellular localization</b>	Nucleus membrane ; Multi-pass membrane protein . Nucleus, nucleolus .
<b>Purification</b>	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
<b>Clonality</b>	Polyclonal
<b>Concentration</b>	1 mg/ml
<b>Observed band</b>	56kD
<b>Human Gene ID</b>	79050
<b>Human Swiss-Prot Number</b>	Q9BVI4
<b>Alternative Names</b>	
<b>Background</b>	similarity:Belongs to the CBF/MAK21 family.,

